


Living up to Life

Leica
MICROSYSTEMS



Leica EM SPF

Self Pressure Freezer



Leica EM SPF: A New Dimension in High Pressure Freezing

The Leica EM SPF Self Pressure Freezer is a system suited for:

- › Cryo-immobilization of biological samples
- › Follow-on procedures for cryo- and room temperature TEM

GENERAL

The Leica EM SPF is based on the principle of Self Pressurized Rapid Freezing introduced by Leunissen and Yi (J. Microsc. 235: 25–35 [2009]). It uses the tendency for water inside a sealed specimen carrier to expand upon cooling, thereby generating pressure intrinsically instead of using an external hydraulic system. This pressure is likely to be the result of crystalline and low density ice formation within the sealed specimen carrier. To achieve pressure (2010 bar) where the melting point of ice is lowered to 251 K (Kanno et al. 1975, Science 189: 880–881 [1975]) 60 % of the water inside the specimen carrier needs to be converted to low density ice.

PRINCIPLE

The Leica EM SPF works with U-tubes as a specimen carrier. The aim is to keep the low density ice located predominantly in the leg areas, whereas the arc area of the U-tube freezes last and while pressurized. This approach for cryo-fixation allows freezing of biological specimens in their native environment without any specific preparation or addition of cryo-protectants, which can alter the initial physiological balance of the sample. Almost any type of cells, free-living bacteria, yeast cells, unicellular organisms etc., can be cryo-immobilized directly after being isolated from their natural habitat.

The Leica EM SPF is a unique entry-level product offering an alternative cryo-fixation method.

YOUR ADVANTAGES

- Enhanced safety for your sample
 - › No cryo-protectant required: Specimens can be cryofixed in their native environment
 - › No Leidenfrost phenomenon: Using propane or ethane as a cryogen prevents the Leidenfrost phenomenon
- Trusted reliability
 - › Control over freezing parameters
 - › Unique U-tube sample carrier
 - › Suitable for samples in suspensions
 - › Defined area of well preserved sample
 - › Reproducible results
- Ease of use
 - › Intuitive operation via touch screen panel
 - › Quick and economical
 - › Compact, mobile and very quiet during freezing





1

5

8

6

7

4

3

2

Leica

Special Features

- 1 Extra large cryo-chamber
providing a comfortable environment for manipulating the U-tubes after cryo-fixation.
- 2 Cryogen container
for ethane or propane
- 3 Platform
for blotting excess propane
- 4 U-tube
(sample carrier)
- 5 Segmenting tool
Designed to make the perfect cut (for precision in segmentation and to reduce mechanical damage during the cutting process)
- 6 Peeling tool
Smooth cutting through the frozen sample and copper with a tungsten carbide knife.
- 7 U-tube separation and temporary storage
- 8 Cryo-transfer container
for CEMOVIS or freeze substitution



Freezing – Operation and Principles

The Leica EM SPF operation is via touch-screen display where every important parameter is adjustable.

Operation

BEND, FILL AND SEAL

U-tubes are made from copper tubes (42 mm long, 0.8 mm OD × 0.4 mm ID) using the bending tool. Filling the U-tube with sample is a standard micropipette procedure after which the open ends are sealed by clamping the leg ends shut using the sealing tool provided.



Principles

UNIQUE U-TUBE CARRIERS

The use of U-tubes in the Leica EM SPF facilitates spatial separation between the area of low density ice formation (predominately in the legs) and the area being pressurized (the arc of the U-tube). The arc of the U-tube is, therefore, the area where the best-preserved specimens are located.

THROUGH THE LOOKING GLASS OF PHYSICS

Low density ice formation causes a volume expansion relative to liquid water. Numerical simulations show that freezing along the tube walls proceeds freezing in the center, producing a strongly curved ice front. This effect is prominent when the immersion speed is the highest. The separation between regions containing low density ice and well-frozen or vitrified parts is considered to be best when the ice front is as flat as possible. This can be achieved by alteration of the freezing parameters for each specific case within the Leica EM SPF program interface. During the immersion movement, the ice formation front inside the tube is always above the level of the cryogen surface. The distance of this ice front depends on the immersion speed.





Details

BENDING THE COPPER TUBE

Bending the copper tube in U-shaped carrier is simply one turn without inducing any stress and deformations in the arc.

OTHER ADVANTAGES


- › LED illumination over the entire work area
- › Thermally controlled injector head
 - ensures conditions for the specimen
 - prevents the arc from pre-cooling
- › Two-step plunging speed






CEMOVIS of *Saccharomyces cerevisiae*

(Courtesy A. Al-Amaudi, DZNE, Bonn)




CEMOVIS of *Pseudomonas deceptionensis*, 30% dextran in PBS


(Courtesy Dr. Carmen Lopez-Iglesias, Lidia Delgado and Elena Mercade, CCI-University Barcelona, Science Park)




Freeze Substitution of *Pseudomonas deceptionensis*

(Courtesy Dr. Carmen Lopez-Iglesias, Lidia Delgado and Elena Mercade, CCI-University Barcelona, Science Park)






Freeze Substitution of *Lingulodinium polyedra*
(*extrusomes*) (Courtesy Elena Lindemann, Fraun-
hofer IGB [Functional Genomics], Stuttgart)



Freeze Substitution of *Lingulodinium polyedra*
(Courtesy PD Dr. Mike Schweikert, University of
Stuttgart)



Freeze Substitution of *Caenorhabditis elegans*
(Courtesy of the Delaware Biotechnology Institute
Bio-Imaging Center (Shannon Modla, Scott Jacobs,
Jeff Caplan and Kirk Czymmek) and University of
Pennsylvania (Jessica Tanis))

WITH THE USER, FOR THE USER

Leica Microsystems operates globally in four divisions, where we rank with the market leaders:

Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

Medical Division

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, "with the user, for the user," describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Groot Bijgaarden	+32	2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1	800 248 0123	847 236 3009
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 5421 2800	3 5421 2896
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People's Rep. of China · Hong Kong	+852	2564 6699	2564 4163
People's Rep. of China · Shanghai	+86	21 6387 6606	21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6779 7823	6773 0628
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Kista	+46	8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44	800 298 2344	1908 246 312
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 236 3009
and representatives in more than 100 countries			