

Living up to Life

Leica
MICROSYSTEMS



Leica EM CPD300

Automated Critical Point Dryer

Controlled Specimen Drying for Biological and Industrial Applications



LEICA EM CPD300 – YOUR BENEFITS

Trusted reliability

- › Fully reproducible processes through full automation
- › Highly reproducible sample preparation

Time saving

- › Reduced process times through new Leica filler/sample holder concept
- › Minimal user interaction time through automation and intuitive software
- › Efficient synergy with Leica pre-treatment instrument EM TP Tissue Processor

Ease of use

- › Ease of use by intuitive software and integrated touch screen
- › Store and recall drying programs
- › Expected process time calculated and displayed according to selected process parameters
- › Timer function
- › Flexibility in sample size
- › Large variety of sample holders



Leica EM CPD300 Critical Point Dryer

The new Leica EM CPD300 dries specimens such as pollen, tissue, plants, insects, etc. as well as industrial samples, for example Micro Electro Mechanical Systems (MEMS), in a fully automated and controlled process. This automated, controlled technique leads to perfect, reproducible results and ensures the same sample quality from every run.

Cost saving

- › Minimized CO₂ consumption through new Leica filler/sample holder concept
- › Minimal user interaction time

Enhanced safety for user

- › Increased safety by software controlled temperature and pressure cut off function
- › Integrated waste separator avoids direct contact

CRITICAL POINT DRYING

The procedure of critical point drying is an efficient method for drying delicate samples for SEM applications. It preserves the surface structure of a specimen which could otherwise be damaged due to surface tension when changing from the liquid to gaseous state. In the past, critical point drying was a time consuming process with low sample reproducibility due to the many manual operations required.

Before drying, many biological samples are commonly prepared through fixation and dehydration and coated after drying with metal such as gold, platinum or palladium to make their surfaces electrically conductive for SEM analysis.

The New Leica Filler Concept – Reduced CO₂ Consumption and Process Time

Sample container with sample holder and filler

The Leica Filler Concept was developed to meet customer demands to reduce CO₂ consumption and process time. With the ability to adjust the chamber volume to the sample size, the exchange volume is reduced, thus minimizing process time and CO₂ consumption.



Filler/holder setting can be programmed to reduce CO₂ consumption and process time.



SAMPLE CONTAINER WITH QUICK RELEASE PIN:

For transferring specimens immersed in exchange fluid to the pressure chamber of the critical point dryer.



HOLDER FOR GRIDS:

Slotted specimen holder for 32 grids. Ideal for drying specimens applied to film-coated grids such as bacteria, viruses, microorganisms, etc.



FILTER DISCS AND POROUS POTS HOLDER:

Specimen holder with four numbered chambers for filter discs, porous pots etc. Ideal for small specimens such as bacteria, viruses, microorganisms, etc. Can also be used without filter discs or porous pots for medium sized samples.



WAFER HOLDER:

Wafer holder for four 2" wafers.



FILLERS:

Filler sizes with a volume of 1/3 and 1/6 of the chamber are available and can be combined with different sample holders and fillers. With the fillers, the chamber volume can be adjusted to the sample size. The exchange volume is reduced which leads to minimized process time and CO₂ consumption.



LEICA EM TP STEM HOLDER:

The EM TP sample basket stem fits directly to the EM CPD300 TP stem holder without unloading the samples from the TP baskets. Can also be used as a universal sample holder for small specimens.



COVERSLIP HOLDERS (can also be used without coverslips for larger samples):

Holder for 22 x 22 mm square coverslips.



COVERSLIP HOLDERS (can also be used without coverslips for larger samples):


Holder for 18 mm diameter coverslips.




COVERSLIP HOLDERS (can also be used without coverslips for larger samples):

Holder for 12 mm diameter coverslips.

Intuitive Software and User Interface with Integrated Touch Screen




Main display before programming: Left side: Temperature, pressure and process time display, including timer function and holder-filler settings; right side: indicates current process status with:



Program display for storing and retrieving programs. Parameter settings for CO₂ IN, fillers, exchange steps and gas OUT settings.



Pop up-filler display: Filler and sample holder settings affect process time. Automated adjustment of the CO₂ IN time, exchange time and gas OUT time to filler-holder status.




Main display: After programming parameters, activated program and process time is indicated (bottom left and top right). Program can be activated with start button.

The Automated Leica SEM Workflow

Leica Microsystems offers a complete preparation package for SEM workflow.

The Leica EM TP is the ideal instrument for tissue processing. After dehydration of the samples with the EM TP, the TP sample basket stem fits directly to the CPD300 TP stem holder without unloading the samples from the TP baskets. Only one user interaction is necessary to facilitate the workflow for dehydration and critical point drying.


Leica Microsystems provides a full range of high quality coating systems as a follow-on procedure to critical point drying.




Comparison of Drying Methods

Water flea (*Daphnia sp.*)


Drying method: **Critical Point Drying with Leica CPD300**



Drying method: **Chemical**






Drying method: **Air**



Courtesy of D. Gruber, University of Vienna,
Austria

Applications


Blood Cells (Erythrocytes and Lymphocytes)




Courtesy of W. Müller, University of Utrecht,
Netherlands

Water flea (*Daphnia* sp.)

(Courtesy of Mag. Daniela Gruber, University of Vienna, Austria)




Tobacco leaf (*Nicotiana tabacum*), trichomes and stomata (Courtesy of Dr. Martin Goldberg and Christine Richardson, University of Durham, U.K.)




Wall Cress (*Arabidopsis thaliana*), stigma

(Courtesy of Dr. Kim Rensing, Applications Specialist Leica Microsystems, University of British Columbia, Canada)






Clawed frog (*Xenopus laevis*), nuclear pores
(Courtesy of Dr. Martin Goldberg and Christine Richardson, University of Durham, U.K.)



Tobacco leaf (*Nicotiana tabacum*), stoma
(Courtesy of Dr. Martin Goldberg and Christine Richardson, University of Durham, U.K.)



Black mold (*Aspergillus niger*), conidiospores
(Courtesy of Dr. Wally Müller, University of Utrecht, The Netherlands)

WITH THE USER, FOR THE USER

Leica Microsystems operates globally in four divisions, where we rank with the market leaders:

Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

Medical Division

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, "with the user, for the user," describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Groot Bijgaarden	+32	2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1	800 248 0123	847 236 3009
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 5421 2800	3 5421 2896
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People's Rep. of China · Hong Kong	+852	2564 6699	2564 4163
People's Rep. of China · Shanghai	+86	21 6387 6606	21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6779 7823	6773 0628
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Kista	+46	8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44	800 298 2344	1908 246 312
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 236 3009
and representatives in more than 100 countries			